OMS-800 series

Operation Microscope

Extensive focus range and high versatility



FTOPCON Healthcare



Bright Wide Field

Topcon's pursuit for perfection is reflected in its continuing development of the OMS-800 range of operating microscopes, adapting them to meet the needs of modern ophthalmic procedures while maintaining the high quality and durablity that made Topcon the world leader in ophthalmic equipment.

Optical Fiber Free Intravitreal Surgery System (OFFISS)

Topcon has developed a state-of-art observation system for vitrectomy procedures that does not require the use of fiberoptic illumination. The Topcon OFFISS lenses avoid complicated focusing by allowing the microscope head and indirect lens to move independently of each other, facilitating a clearly focused image at all times. The image inverter activates automatically whenever the OFFISS is in use. The indirect lens can quickly and simply be exchanged for another, saving time and increasing efficiency.



Small 40D Lens¹ (Outer Diameter φ23) (Outer Diameter φ28) (Outer Diameter φ19.4) (Outer Diameter φ24)







120D Lens



Small 120D Lens² (Outer Diameter \$18)



Anterior Lens³ (Outer Diameter $\phi34$)

Small diameter 40D lens Small diameter 120D lens Anterior eye section observation lens

TOPCON OFFISS LENSES

Small 40D Lens

- · Assists with membrane peeling in the macular region.
- · Beneficial for highly myopic eyes by avoiding contact between the surgical tools and lens.
- · Provides a crisp, wide-angle view with remarkable stereopsis, giving a clear view of the posterior pole.4

40D Lens

- · Bright, stereoscopic view.
- · Ideal for posterior segment procedures.
- · Combined with the microscope illumination, the characteristics of the lens make the use of additional fiber optic illumination unnecessary, thereby enabling

bimanual procedures and saving time.

80D Lens


- · Allows observation from the posterior segment out to the intermediate peripheral zone of the retina.
- · Can be used in combination with fiber optic illumination.

120D Lens

- · Useful for vitreous surgery and photocoagulation of the central and peripheral areas up to the Ora Serrata.
- · Provides a field of view of 130 degrees with good stereopsis.
- · Can be used under air substitution in combination with a wide-angle fiber optic endoilluminator.


Small 120D Lens


- · Compact lens takes up minimal space in the operating field.
- · Does not interfere with the use of surgical instruments.
- · Wide angle view of 100 degrees can expand up to approximately 130 degrees with the use of air substitution.




Small 40D Lens allow better vertical access for surgical tools than the conventional 40D lens.









FOR CATARACT AND VITREOUS SURGERY

Superb image quality for cataract and vitreous surgery

With the advancement of cataract surgery and phacoemulsification techniques, an increasing number of surgeons are performing simultaneous cataract and vitreous surgeries. By using a three-mode illumination system, the OMS-800 provides an improved red reflex with better shadow and contrast, even under conditions of low illumination.



Continuous Curvilinear Capsulorhexis (CCC)



Hydrodissection




Phacoemulsification and Aspiration (PEA)


Easy switch among illumination modes

Three different illumination modes are available to meet all surgical lighting needs.


Different modes are easily selected using the footswitch.



In this mode, the illumination, brightness, stereoscopic view, and shadow contrast are perfectly balanced for superior observation clarity. The illumination is always optimum regardless of the position of the patient's eye.



This illumination mode generates a particularly good red-reflex and is very useful during anterior capsulotomies.



The combination of illumination and yellow filter is particularly advantageous during long procedures to prevent phototoxicity.

ENHANCED FUNCTIONALITY

Low intensity illumination enables clear observation while preventing light damage

Superbly designed optics provide optimum illumination, eliminating harmful wavelengths and unnecessary brightness. Low light intensity also helps to prevent light-related damage to the retinal tissues.

Coarse focusing*

The coarse focusing mechanism allows the optical head to be quickly elevated during surgery and then brought back to the desired working position. This feature is particularly useful during IOL insertion and other procedures that momentarily require more space between the patient and the microscope.

* Refer to the component list for further detail.

Apochromatic optics

The optics of the system are designed to greatly limit the effect of chromatic aberration.

Anti-stain coating

The OMS-800 employs an anti-stain coating - the optical components remain clear and maintain their quality for a longer period.

Multi-function footswitch enhances operating efficency

The multi-function footswitch permits surgeons to control virtually all of the OMS-800 functions without removing their hands from the operative field. Without any hand movement, the surgeon can adjust the illumination, zoom magnification, focus, illumination angle, and X-Y positioning.

The control layout on the footswitch is conveniently arranged in the most accepted configuration.

Comfortable operating posture

The ergonomically designed optical head with a built-in beam splitter and adjustable eye pieces allows the surgeon to maintain a comfortable posture throughout the surgery. A key component to this comfort is the variable angle binocular tubes that allow for the setting of a personal viewing position from 45 to 90 degrees. This flexibility ensures a comfortable operating stance even when using OFFISS.

Comfortable working distance

The OFFISS system provides an comfortable working distance between the OFFISS lenses and the patient's eye.

Easy bulb exchange

The illumination bulb can be simply replaced by accessing the light housing using a rotating lever.

A warning lamp indicates when the spare lamp is burned out to ensure there is always an operational bulb available.

Electromagnetic locking system

The optical head can be quickly and accurately positioned for surgery and held in place by the fast-acting electromagnetic locking system. (OFFISS/Pro)



OPTIONAL ACCESSORIES

OFFISS lens set

Standard components include: front lens holder, anterior segment observation lenses, 40D, small 40D, 80D, 120D, and small 120D. The boxes and lenses are easily maintained using an autoclave.



Assistant microscopes

The assistant microscope provides an additional viewer with bright, crisp images on the same visual axis as those seen by the main surgeon.

Coaxial binocular tube

This tube allows observation of the same field of view and magnification as the main surgeon. The angle of the binocular eyepieces is adjustable from 30 to 110 degrees, offering the assisting surgeon a comfortable viewing angle.



Coaxial binocular tube

O degree assistant microscope

The angle of the binocular eyepieces is adjustable from 45 to 90 degrees, offering the assisting surgeon a comfortable viewing angle. In addition, a separate focus adjustment is available for the assistant surgeon.



O degree assistant microscope

TV relay lens

The compact TV relay lens permits the attachment of a CCD camera useful for documentation and teaching. The relay lens accepts the most popular 1/2" and 1/3" CCD cameras with C mount or bayonet mount and is easily connected to the OMS-800.



The image is as sample. Please check with your local distributor in details.

OMS-800 OFFISS

OFFISS offers a scope of possibilities for vitreoretinal surgery. Equipped with the OFFISS lenses mechanism, electromagnetic brakes, and sophisticated electronics, this model is the highest specification for intravitreal surgery, as well as other ophthalmic procedures.



Electromagnetic brakes and sophisticated electronics confer the OMS-800 Pro the flexibility to facilitate virtually any type of ophthalmic surgical procedure.

OMS-800 Standard

Equipped with most of the state-of-the-art features of the OMS-800 range, the OMS-800 Standard answers the need for a simpler, user-friendly operation microscope. Manual brakes and ease of mobility make the OMS-800 an affordable yet advanced unit for all ophthalmic procedures.





Component

	OFFICE	Due	Chandand
	OFFISS	Pro	Standard
OFFISS	V	-	-
Electromagnetic locking	V	V	-
Coarse focusing	V	V	-
Inverter	V	-	-
Apochromatic optics	V	V	V
Beam splitter	V	V	V

Specifications

Microscope type	Galileo type	
Magnification change type	Electric zoom continous change	
Eyepiece (Eyepiece magnification)	12.5x	
Objective lens	f=200mm	
Display magnification(x)	4.2/5/6/7/8/9/10/11/13/15/17/19/21	
Total magnification	4.2x - 21x	
1st arm length (Distance between shafts)	375mm	
1st arm rotation range	300°	
2nd arm length (Distance between shafts)	OFFISS, Pro : 990mm Standard : 875mm	
2nd arm rotation range	300°	
2nd arm vertical movement range	600mm	
2nd arm mounting weight	OFFISS, Pro : 6kg - 18kg Standard : 9kg - 21kg	
Power supply	AC 100-240V, 50-60Hz	
Power supply input	280VA	
Dimensions		
Base (Base unit)	720mm(W) x 720mm(D)	
Base (Base total height)	1,865mm	
Weight	OFFISS : 250kg Pro : 247kg Standard : 244kg	
Permitted weight for accessories	OFFISS : 4.8kg Pro : 6.8kg Standard : 7.3kg	

Photos Courtesy of: Professor Masayuki Horiguchi, MD Ophthalmology Department Fujita Health University

Associate Professor Kiyoshi Suzuma, MD Department of Ophthalmology and Visual Science Graduate School Biomedical Science, Nagasaki University

IMPORTANT
In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation. Not available for sale in all countries. Please check with your local distributor for availability in your country.
Not all products, services or offers are approved or offered in every market, and products vary from one country to another.
Contact your local distributor for country-specific information and availability.







75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN. Phone: +81-(0)3-3558-2522/2502 Fax: +81-(0)3-3965-6898 https://topconhealthcare.jp

